Test Report

Perceptra Inspectra CXR Version 3

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
Inspectra CXR is an AI-based software that analyzes adult frontal chest x-rays for signs of abnormalities in order to prioritize them for clinical review. The device aims to assist in prioritizing and triaging chest x-ray images. Inspectra CXR identifies features suggestive of 8 abnormality conditions:

1. Pleural effusion
2. Tuberculosis
3. Lung opacity
4. Pulmonary edema
5. Mass
6. Nodule
7. Cardiomegaly
8. Atelectasis

The standalone software consists of both on-premise and on-cloud modules. The on-premise module interacts with the hospital Picture Archiving and Communication System (PACS) in order to retrieve input DICOM images and de-identify data before sending it to the on-cloud module to process. The on-cloud deep learning algorithms then analyze chest x-rays and output structured reports detailing probabilities of abnormalities in the image. The software is capable of producing a secondary capture DICOM image with AI analysis results to send back to the PACS system or display on MyInspectra web application. The software flags cases with suspected findings, showing passive notification on the work list to notify qualified physicians to urgently review the suspected cases. X-rays without an identified anomaly are placed on the worklist for a routine review, which is the current standard of care. Inspectra CXR does not mark, highlight, or direct user attention to a specific location on the original chest x-ray image. The device generates a secondary image from the original image with an overlay that has to be turned on by qualified physicians when additional information is needed in the clinical review.

The system is compatible with adult chest radiographs starting from 15 years old. The system sorts out pediatric chest radiographs, Posterior-Anterior (PA) Upright radiographs, irrelevant radiographs, and non-digital chest radiographs before processing. The system provides preliminary screening which allows only images indicated as the lung and with the patient's age stated. Images without the part of the patient's body identified or without the age stated will not be processed. The device intends to only assist radiologists /trained clinicians to facilitate clinical review of chest x-rays. It does not include clinical or patient historical data or recommends treatments or diagnoses. The final decision has to be provided by attending physicians or radiologists.

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ICC/Kappa Statistic	Strength of Agreement
<0.00	Poor
0.00 – 0.20	Slight
0.21 – 0.40	Fair
0.41 – 0.60	Moderate
0.61 – 0.80	Substantial
0.81 – 1.00	Almost Perfect

¹Landis, J. R., & Koch, G. G. (1977). The Measurement of Observer Agreement for Categorical Data. In *Biometrics* (Vol. 33, Issue 1, p. 159). JSTOR. <https://doi.org/10.2307/2529310>

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Finding	N	Threshold	Pairwise Agreement		Cohen's Kappa	
			"B" vs "B"	"B" vs AI	"B" vs "B"	"B" vs AI
Abnormality	501	*	0.9007	0.8811	0.7993	0.7578
Tuberculosis	422	0.43	0.9635	0.9411	0.9261	0.8817
Opacity [#]	464	0.38	0.9087	0.9244	0.8190	0.8488
Pleural Effusion	99	0.17	0.9578	0.9544	0.7838	0.7704

[#]Opacity combines small and large opacities.

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
Finding	Pairwise Agreement			Cohen's Kappa		
	B vs "B"	"B" vs AI	Diff	"B" vs "B"	"B" vs AI	Diff
Abnormality	0.9007	0.8811	-1.96%	0.7993	0.7578	-4.15%
Tuberculosis	0.9635	0.9411	-2.24%	0.9261	0.8817	-4.44%
Opacity [#]	0.9087	0.9244	+1.57%	0.8190	0.8488	+2.98%
Pleural Effusion	0.9578	0.9544	-0.34%	0.7838	0.7704	-1.34%

[#]Opacity combines small and large opacities.

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	Minimal Requirements		Optimal Requirements	
	Sensitivity	Specificity	Sensitivity	Specificity
Smear-replacement test	Overall >80%	98%	Overall >95%	98%
	Positive >99%		Positive >99%	
	Negative >60%		Negative >68%	
Non-sputum based biomarker test	Overall >65%	98%	Positive >98%	98%
	Positive >98%		Negative >68%	
Triage test	90%	70%	95%	80%

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